Sampling & Testing

QuarkE-MIC-PM

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Sampling & Testing

QuarkE-MIC-PM

Planktonic Microorganisms Sampling and Testing Kit


The QuarkE-MIC-PM system is based on the Most Probable Number (MPN) method and is dedicated to the rapid sampling and testing of MIC-prone microorganisms such as SRB, TGB, and FB in fluid media within equipment and facilities in marine and terrestrial oil and gas fields. This system and its user manual refer to standards such as NACE TM0194 and NACE TM0212.

User Manual:


Step1: When using this test kit on-site, safety helmets, anti-static work clothes, work boots, safety glasses, earplugs, gloves, and hydrogen sulf respirators should be worn, or configured according to the on-site safety management regulations.

Step2: Select representative sampling points, for facilities such as oil and water pipelines, the sampling points should be located at the 6 o'clock position on the horizontal section. Use the disinfectant cloth in the test kit to thoroughly disinfect the exterior of the sampling point.

Step3: Open the sampling point valve and displace the sample at the sampling point for 1-3 minutes, then close the valve. Open the sterile sampling bag in the test kit, connect the bag opening to the sampling point, quickly open the sampling point valve slightly, and collect 200 ml of sample at a low flow rate, and quickly tighten the cap.

Sampling & Testing


QuarkE-MIC-PM

User Manual:


Step4: Record the sampling date, time, pressure, temperature, pH (if no data is available, it needs to be tested on-site), dissolved oxygen, 2S concentration, other water chemistry data that can be collected, sample color (especially suspended metal sulfides or black water), turbidity, odor, presence of or precipitate, flow rate, and the daily production of oil, gas, and water.


Step5: Place the sterile bag containing the collected sample at rest (resting time should not exceed 1h). After the oil and water have separated, use syringe provided in the test kit to draw 1ml of the water sample, and quickly inject it into the No. 1-1 culture bottle. Shake liquid in the culture bottle evenly.

Step6: Replace the syringe and slowly and evenly withdraw 1 ml of liquid from 1-1 and inject it into 1-2. Repeat this process, the solution is diluted to the specified multiple (for the planktonic state, it is recommended to dilute 2*3, for the attached state and the state within the corrosion products, it is recommended to dilute 2*6).

Step7: Record the sampling location, sampling time, and inoculation time on the culture medium packaging box, and place the culture medium in a constant temperature box (temperature of the constant temperature box should be set according to the temperature of the sampled equipment). Take out and read the results after 7 days and 14 days, and take positive photos for record.

Product List:

1	Gloves	1
2	Alcohol Disinfectant Tablets	2
3	Sterile Sampling Bag	2
4	Sterile Syringe	12
5	High Precision Culture Medium	12

Contact Us:

TEL: 0775-23228012

MP: 134 3786 2586

E-mail: info@qeelab.com

Add.: Shenzhen, Guangdong